CLINIPET⁺

Autoclavable

ENGLISH 1 – 14

DEUTSCH 15 – 28

FRANÇAIS 29 – 42

ESPAÑOL 43 – 56

POLSKI 57 – 70

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


CE IVD

HTL
HTL LAB SOLUTIONS

PZ HTL S.A. is in possession of a registered and certified quality management system which includes the development, production, and sales of high quality Liquid Handling products.

1 - INTRODUCTION

CLINIPET+ pipettes are designed for accurate and repeatable measure and transfer of liquids of constant volume from 2 to 1000 μ l.

Liquid is aspirated into a disposable tip attached to the pipette shaft.

NOTE: Disposable tips ensure maximum safety and eliminate possibility of cross-contamination of liquid samples.

CLINIPET+ pipettes have a blow-out mechanism ensuring complete removal of liquid from the tip.

CLINIPET+ pipettes have a built-in tip ejector. Using it makes removing the tip from the pipette shaft easier, and protects against contact with a contaminated tip. When using narrow tubes it may be necessary to remove the tip ejector. It is simply removed by pulling down.

The accuracy (A) and repeatability (P) of liquid volume depends on the quality of tips used. The values for accuracy and repeatability were obtained using the **HTL** tips. Only the use of **HTL** tips guarantees accuracy and repeatability of liquid sampling.

Figure 1 shows the outside elements of the pipette and their materials:

A. Pipetting pushbutton	Polypropylene (PP)
B. Eject pushbutton	Polypropylene (PP)
C. Handle	Polypropylene (PP)
D. Ejector	Polypropylene (PP)
E. Pipette shaft	Polyvinylidene fluoride (PVDF)

Please, get acquainted with the figure and its description before reading the instruction manual.

CONTENTS

1 - INTRODUCTION
2 - ASPIRATION AND DISPENSING LIQUIDS
3 - HTL PIPETTING TIPS
4 - RECOMMENDATIONS
5 - RECALIBRATION
6 - PIPETTING DENSE AND VISCOUS LIQUID
7 - TROUBLESHOOTING
8 - CLEANING AND STERILISATION
9 - PIPETTE KIT
10 - SPARE PARTS

Pipette model	Cat. number	Volume [μ l]	Accuracy A [%]	Precision P [%]	Volume change ΔV [μ l]	Tip volume	Shaft colour
CP 2	2202	2	± 3.0	± 1.5	0.7	200 μ l	yellow
CP 5	2204	5	± 2.0	± 1.0			
CP 7	2206	7	± 1.5	± 0.7			
CP 10	2208	10	± 1.2	± 0.5	1.3		
CP 15	2210	15	± 1.0	± 0.5			
CP 20	2212	20	± 0.9	± 0.4	1.7		
CP 25	2214	25	± 0.9	± 0.4			
CP 30	2216	30	± 0.9	± 0.3	3.5		
CP 40	2218	40	± 0.9	± 0.3			
CP 44.7	2220	44.7	± 0.9	± 0.3			
CP 50	2222	50	± 0.9	± 0.3	7.0		
CP 60	2224	60	± 0.9	± 0.3			
CP 70	2226	70	± 0.9	± 0.3			
CP 75	2228	75	± 0.9	± 0.3	13.0		
CP 80	2230	80	± 0.9	± 0.3			
CP 90	2232	90	± 0.8	± 0.3			
CP 100	2234	100	± 0.8	± 0.3	31.0		
CP 120	2236	120	± 0.8	± 0.3			
CP 150	2238	150	± 0.7	± 0.3			
CP 200	2240	200	± 0.6	± 0.2	60.0		
CP 200A	2242	200	± 0.6	± 0.2			
CP 220	2244	220	± 0.6	± 0.2			
CP 250	2246	250	± 0.6	± 0.2	1000 μ l	blue	
CP 300	2248	300	± 0.6	± 0.2			
CP 400	2250	400	± 0.6	± 0.2			
CP 450	2252	450	± 0.6	± 0.2			
CP 500	2254	500	± 0.6	± 0.2			
CP 600	2256	600	± 0.6	± 0.2			
CP 700	2258	700	± 0.6	± 0.2			
CP 750	2260	750	± 0.6	± 0.2			
CP 800	2262	800	± 0.6	± 0.2			
CP 900	2264	900	± 0.6	± 0.2			
CP 1000	2266	1000	± 0.6	± 0.2			

The accuracy and precision figures are obtained using HTL tips, using a gravimetric method, performing at least 10 measurements of distilled water at the temperature of $20 \pm 1^\circ\text{C}$, according to EN ISO 8655 standard.

The pipette design enables the user to perform the recalibration process according to the rules presented in section 5.

Volume change ΔV is the change of the distilled water volume aspirated by the pipette caused by a single complete turn of the calibration key - 24 increments.

2 - ASPIRATION AND DISPENSING LIQUIDS

Attaching the tip

- Select an appropriate tip (see Section 3).
- Press the ejector pushbutton with your thumb from the top, (Fig. 2A).
- Pressing the ejector pushbutton, push the pipette cone into tip opening, (Fig. 2B).

NOTE! Never aspirate liquid into a pipette without a tip attached.

Pre-rinsing the tip

When a tip will be used for the first time, especially when pipetting liquids of viscosity or surface tension different from water (e.g. serum, organic solvent), inner surface of the tip is wetted with the liquid remaining inside. Wetting inner surfaces causes errors, as the volume of the first dispensed sample will be smaller than that of successive samples. To avoid those errors, the tip should be pre-rinsed with the pipetted liquid before the first sample is dispensed. This is accomplished by aspirating a sample of the required volume and dispensing it back. When the tip is pre-rinsed, the next pipetted volume will be free of any wetting error.

Aspiration of liquid

- Press the pushbutton slowly and smoothly to the first stop, (Fig. 3A).
- Holding the pipette vertically, immerse the tip into the sampled liquid. The tip immersion depth depends on the tip used, and the sampled volume.

Tip [μ l]	Immersion depth [mm]
200	2÷3
1000	2÷4

- Release the pipetting pushbutton slowly and uniformly, aspirating liquid into the tip, (Fig. 3C).
- Wait one second and withdraw the tip from the liquid.

Liquid Dispensing

- Hold the pipette slightly inclined ($10 \div 40^\circ$ off the vertical), and place the tip orifice against the inside wall of the receiving vessel.

- Press the pipetting pushbutton smoothly and slowly to the blow-out stop, dispensing the liquid, (Fig. 3D).
- Wait one second.
- Press the pipetting pushbutton to the final stop, blowing out any remaining liquid.
- Hold the pipetting pushbutton depressed, and remove the tip from the receiving vessel, drawing it against the vessel inner wall.
- Release the pipetting pushbutton to its starting position, (Fig. 3F).

Replace the tip with a new one whenever a different liquid is to be pipetted.

Ejecting the Tip

- Place your thumb against the ejector pushbutton side (Fig. 2D), press and turn it. The tip will be ejected from the pipette shaft.

3 - HTL PIPETTING TIPS

HTL tips are made of pure polypropylene. Their quality ensures proper operation of **CLINIPET+** pipettes, and secures pipetting accuracy and repeatability.

Accuracy (A) and repeatability (P) of **CLINIPET+** pipettes was determined using **HTL** tips. Using other tips may result in deterioration of pipetting accuracy and repeatability.

Tip 200

This tip is used for volumes from 2 μl to 200 μl .

Tip 1000

This tip is used for volumes from 200 μl to 1000 μl .

4 - RECOMMENDATIONS

Observing the following recommendations will ensure maximum possible accuracy and precision of liquid sampling.

- Operate the pipette pushbutton slowly and smoothly when liquid aspirating and dispensing.
- The depth of tip immersion in the aspired liquid should be kept at necessary minimum, and should remain constant during aspiration.

- While operating, the pipette should be held in a vertical position.
- The tip should be replaced with a new one each time a different kind of liquid is handled.
- The tip should also be changed if visible droplets of liquid remain inside it.
- Each new tip should be pre-rinsed with the pipetted liquid.
- Pipetted liquid must not enter the pipette shaft. To ensure this:
 - operate the pipetting pushbutton slowly and smoothly;
 - do not lay the pipette down in a horizontal position if there is any liquid in the tip.
- Do not aspirate liquids of temperatures above 70°C.
- After pipetting acids and aggressive liquids, it is recommended to disassemble the pipette and rinse the pipette plunger, shaft, seal and other elements in distilled water.

5 - RECALIBRATION

HTL pipettes are calibrated by gravimetric method, using **HTL** tips and distilled water, at temperature $20 \pm 1^\circ\text{C}$ according to EN ISO 8655 standard.

If during pipette operation or checking you find that the accuracy error (A) - the difference between the real aspirated volume and the nominal pipette volume, exceeds the permissible value given in the table in section 1, the pipette recalibration procedure should be carried out.

Before starting the recalibration procedure check the conditions under which the accuracy error (A) has been determined and ensure whether:

- the ambient temperature and the temperature of the pipette, tips and water is the same
- the density of the liquid used is close to that of distilled water
- balances of appropriate sensitivity have been used

Volume checked [μ l]	Balance sensitivity [mg]
2 - 10	≤ 0.001
10 - 100	≤ 0.01
> 100	≤ 0.1

- the conversion factor mg/ μ l has been taken into consideration in the calculations
- the requirements given in sections 2 and 4 have been fulfilled.

If you follow the instructions and the accuracy error (A) still exceeds the permissible value carry out the pipette recalibration procedure again. Carry out recalibration procedure in the following way:

- use distilled water for measurements
- obey the rule that ambient temperature and the temperature of the pipette, tips and water should be the same within the range 20-25°C and stabilised during weighing within $\pm 0.5^\circ\text{C}$
- perform 5 aspirations, weigh each time using the balance which meets the sensitivity requirements
- calculate the average value of the aspirations
- calculate the average aspirated volume, in μ l, multiplying the average aspiration amount in [mg] by the distilled water density coefficient, [mg/ μ l] which depends on temperature and pressure according to the following table:

Density coefficient values for distilled water

Temperature [$^\circ\text{C}$]	Pressure [kPa]		
	95.0	101.3	105.0
20	1.0028	1.0029	1.0029
21	1.0030	1.0031	1.0031
22	1.0032	1.0033	1.0033
23	1.0034	1.0035	1.0036
24	1.0037	1.0038	1.0038
25	1.0039	1.0040	1.0040

- determine the difference between the average aspirated volume and the nominal volume of the pipette.

If the difference does not exceed the permissible value of the accuracy error (A), the pipette is set properly and does not require recalibration.

Otherwise:

- remove the pipette pushbutton, (Fig. 4A)
- insert the calibration key (with the smaller diameter end) into the pipette opening, the key should enter the calibration screw located inside the pipette, (Fig. 4B)
- to reduce the aspiration volume turn the key clockwise, to increase the aspiration volume turn the key counter-clockwise, (Fig. 4C). One full turn of the key corresponds to the change of the aspirated volume of distilled water by the amount ΔV given in the table in section 1 - 24 increments
- take the calibration key out, fix the pipette pushbutton and determine the average aspiration amount again, (Fig. 1D)
- the determined volume should be within the permissible range accuracy error (A)
- otherwise, repeat the recalibration procedure once more.

More information on the recalibration procedure can be found on www.htl.com.pl.

6 - PIPETTING DENSE AND VISCOUS LIQUID

The specification of accuracy and repeatability of pipetting with **CLINIPET⁺** pipettes is based on sampling distilled water. When pipetting other liquids of different physical qualities like density, viscosity and surface tension, it is recommended to aspirate and dispense that liquid much slower.

When aspirating, the tip should be kept in the pipetted liquid for at least 2 seconds. Similarly, when liquid dispensing, you should wait a minimum 2 seconds before you press the pipetting pushbutton to its ultimate end.

If the difference between the real volume of aspirated liquid and the nominal volume of the pipette is considerable and inadmissible, the pipette recalibration procedure can be carried out for the liquid actually aspirated.

Carry out the recalibration procedure as for distilled water by changing the liquid only.

For calculating the average volume, in l, aspirated by the pipette, use the following formula:

$$\text{Average aspirated volume } [\mu\text{l}] = \frac{\text{Average aspirated amount } [\text{mg}]}{\text{Liquid density } [\text{mg}/\mu\text{l}]}$$

Mark the liquid the pipette recalibration was carried out for.

Remember that in case of aspirating a liquid of a density close to that of distilled water you should recalibrate the pipette again.

7 - TROUBLESHOOTING

If you notice an improper pipette operation, find the reason and eliminate the fault. Doing this, follow the instructions in the given sequence. Replacement of elements into new ones may happen only exceptionally, and should not occur under normal pipette use.

Droplets of liquid remain in the pipette tip.

- The pipette is emptied too fast.
Decrease the speed of pressing the pipette pushbutton.
- The tip wettability has increased due to extensive use.
Replace the tip with a new one.

Droplets of air appear in the aspirated liquid.

- The pipette tip immersion is too shallow.
Immerse the tip deeper according to the instructions.
- Liquid is aspirated too fast.
Decrease the speed of aspiration.
- The pipette tip is inadequately pressed onto the pipette cone.
Press the pipette tip harder.
- The tip is damaged or worn out due to extensive use.
Replace the tip with a new one.

Liquid drops out from the tip.

- The pipette tip is inadequately pressed onto the pipette cone.
Press the pipette tip harder.

- The pressure plate (10) of the nut (14) is loose, (Fig. 6).
Remove the ejector (16), pressing the ejector pushbutton (2) unscrew the nut (14), pull out the plunger (8) and plunger spring (9). Tighten the pressure plate (10).
- The pipette inner parts are dirty because of prolonged pipetting of corrosive agents or because liquid got inside the pipette.
*Remove the ejector (16), pressing the ejector push-button (2) unscrew the nut (14), pull out the plunger (8) and plunger spring (9). Unscrew the pressure plate (10) from the nut (14), and remove all remaining elements from the nut (14) and the shaft (15). Wash the shaft, nut, plunger and all other elements in distilled water or isopropyl alcohol, and dry them. Grease the gasket (12) and O-ring (13) with lubricant that has been included with each pipette.
*Re-assemble the pipette in the reverse order.**
- The shaft (15) sealing surface is damaged.
*Disassemble the pipette as described above. Wash and dry the elements. Grease the gasket (12) and O-ring (13) with lubricant.
Replace the shaft (15) with a new one.
Re-assemble the pipette in the reverse order.*
- The gasket (12) and/or O-ring (13) is damaged.
*Disassemble the pipette as described above. Wash and dry the elements. Replace the gasket (12) and O-ring (13) and grease them with lubricant.
Re-assemble the pipette in the reverse order.*
- The plunger (8) is damaged.
*Disassemble the pipette as described above. Wash and dry the elements. Grease the gasket (12) and O-ring (13) with lubricant.
Replace the plunger (8) with a new one.
Re-assemble the pipette in the reverse order.*
Replacement of the plunger requires recalibration of the pipette.

Exceeded value of the accuracy error (A)

If you find that the accuracy error (A) - the difference between the real volume aspirated and the

nominal volume of the pipette, exceeds the permissible value given in the table in section 1, you should recalibrate the pipette according to the rules presented in section 5.

Increased ejector pushbutton turning resistance.

- The handle (5) inner surface is dirty.

Remove the ejector (16), pressing the ejector push-button (2) unscrew the nut (14), pull out the ejector spring (7). Slide the body (6) out from the handle (5).

Clean the handle inner surface and the body outer surface with dry cotton wool. Grease the ejector spring, handle inner surface and the body outer surface which enters the handle, with lubricant.

Re-assemble the pipette in the reverse order.

Attention: Use only lubricant that has been included with each pipette.

If the pipette malfunction persists, deliver it to your **HTL** service representative.

Before returning the pipette make sure that it is free from any chemicals, radioactive or microbiological contamination that might pose any danger during transportation and servicing.

If possible, clean the pipette to perform the sterilisation.

8 - CLEANING AND STERILISATION

Cleaning:

External surfaces of the pipetting pushbutton (1) - Fig. 6, the ejector pushbutton (2), the handgrip (5) and the body (6) can be cleaned with tissue dipped in isopropyl alcohol. Remaining parts removed from the pipette during disassembly (7 to 17) can be washed with distilled water or isopropyl alcohol.

Sterilisation:

The pipette can be subjected to sterilisation, as a whole, in the autoclave at temperature 121°C for 20 minutes. After performing the sterilisation the pipette should be dried and cooled down to room temperature.

It is recommended:

- to sterilise the pipettes in an autoclave with initial vacuum and drying cycle,
- to test pipette calibration every 10 sterilisation cycles.

The pipettes are UV resistant which was confirmed by our tests. The recommended distance from the radiation source to exposed element should be not less than 50 cm. Prolonged and very intense UV exposure can cause de-colouration of pipette parts, without affecting its performance.

9 - PIPETTE KIT

Pipettes are delivered in the kits including:

- Pipette,
- Instruction manual,
- Calibration key,
- Pipette stand,
- Lubricant.

The way of fixing the stand is shown in Fig. 5.

10 - SPARE PARTS

Fig. 4, 6 show the pipette design and its elements.

Replacement parts are marked with bold type-faces.

- 1. Pipetting pushbutton**
2. Ejector pushbutton
3. Body cover
4. Handle cover
5. Handle
6. Body
- 7. Ejector spring**
- 8. Plunger**
- 9. Plunger spring**
- 10. Pressure plate**
- 11. Pressure sleeve**
- 12. Gasket**
- 13. O-ring**
- 14. Nut**
- 15. Shaft**
- 16. Ejector**
- 17. Pressure sleeve spring**
- 18. Calibration key**

All the spare parts indicated in Fig. 4, 6 can be ordered from your **HTL** representative (type of pipette and name of the part for this pipette should be specified).

Warning: The replacement of the plunger requires conducting of calibration procedure according to section 5.

All rights reserved. Product described in this manual is subject to availability and technical modification. Errors excepted. PZ HTL S.A. reserves the right to improve, enhance or otherwise modify its products without prior notification.

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© 2015 PZ HTL S.A.