



BioMedical

CryoBiological Freezer and Dewar Accessories Catalog



Table of Contents

Transfer of Liquid Nitrogen

| | |
|-------------------------------|---|
| pouring spouts..... | 3 |
| manual discharge devices..... | 3 |
| swivel dipper..... | 3 |
| dipper..... | 4 |
| phase separator..... | 4 |
| transfer hoses..... | 4 |
| union coupler..... | 5 |

Shipping Accessories

| | |
|----------------------------|---|
| cryosleve..... | 5 |
| cabosil blocks..... | 5 |
| square cabosil insert..... | 6 |

IATA Accessories

| | |
|--|---|
| secondary container..... | 6 |
| screws for secondary container..... | 6 |
| absorbant paper for secondary container..... | 7 |
| cryovials..... | 7 |

Protective Shipping Cartons

| | |
|--------------------------|-----|
| shipping containers..... | 7-8 |
|--------------------------|-----|

Roller Bases

| | |
|---------------------------------------|---|
| roller bases for aluminum dewars..... | 8 |
|---------------------------------------|---|

Corks & Covers

| | |
|----------------------------------|---|
| corks & covers for freezers..... | 8 |
|----------------------------------|---|

Supply Accessories

| | |
|-------------------------------|---|
| AI measuring stick..... | 9 |
| freezer measuring sticks..... | 9 |

Storage in Canes & Canisters

| | |
|----------------------------------|----|
| straw canisters..... | 9 |
| center canister..... | 10 |
| C-1 cane..... | 10 |
| C-2 cane..... | 10 |
| C-7 cane..... | 11 |
| divider and center canister..... | 11 |
| cabosil canister insert..... | 11 |

Storage in Boxes

| | |
|---------------------------------|----|
| plastic mini boxes..... | 12 |
| large plastic vial boxes..... | 12 |
| large cardboard vial boxes..... | 12 |

Storage in Racks

| | |
|---------------------------|----|
| square canister rack..... | 13 |
| SUC box..... | 13 |
| SUC marker cane..... | 13 |
| SUC cane marker..... | 14 |

Electronic Alarms

| | |
|--------------------------------|----|
| TEC 3000..... | 14 |
| MVE datalogger..... | 14 |
| Planer datalogger..... | 15 |
| Bat-1-B..... | 15 |
| liquid level alarm (240V)..... | 15 |

Freezer Accessories

| | |
|-----------------------------------|----|
| 1st generation converter kit..... | 16 |
| 2nd generation converter kit..... | 16 |
| tank switcher..... | 16 |
| tee valve assembly..... | 17 |
| vapor platform..... | 17 |
| power conditioner..... | 17 |
| freezer step assembly..... | 18 |
| temperature sleeves..... | 18 |

Safety

| | |
|----------------------------|----|
| apron..... | 18 |
| cryo-gloves..... | 19 |
| cool reach / cryoclav..... | 19 |
| epoxy kit..... | 19 |

Freezer Inventory Systems

CryoBio Accessory Catalog

Pouring Spout



| | | |
|---------|--------------|---------------|
| 9720729 | 2.2" Neck ID | Lab 5, 10, 20 |
| 9711909 | 2.0" Neck ID | old Lab 20 |
| 9711929 | 2.5" Neck ID | Lab 30, 50 |

Manual Discharge Device



| | | |
|-------------------|--------|------------|
| 10668101 | 2.2"ID | Lab 10 |
| 13467208 | 2"ID | old Lab 10 |
| 13464867 | 2"ID | old Lab 20 |
| (prior to 3/2007) | | |
| 13484139 | 2.2" | new Lab 20 |
| (after 3/2007) | | |
| 10668004 | 2.5"ID | Lab 30 |
| 10668021 | 2.5"ID | Lab 50 |

Swivel (Picasso) Dipper



| | | |
|---------|----------|-----------|
| 9711679 | 21.25"lg | 17ml cap. |
|---------|----------|-----------|

CryoBio Accessory Catalog

Dipper



45ml approx 1.5 ounces

| | | |
|---------|-----------|------------|
| 9711569 | 16.25"lg | Lab 5 |
| 9711589 | 19.625"lg | Lab 10 |
| 9711619 | 21.75"lg | Lab 20, 30 |
| 9711669 | 28.25"lg | Lab 50 |

103ml approx 3.5 ounces

| | |
|----------|----------|
| 11555669 | 22.75"lg |
|----------|----------|

Phase Separator



| | | |
|----------|--------|----|
| 10615885 | small | Al |
| 10615877 | medium | Al |
| 10615869 | large | Al |

Transfer Hose



LN2 transfer hose 1/2" ODT

| | |
|---------|---------|
| 9713159 | 4' long |
| 9713109 | 6' long |

Combos

| | |
|----------|-------------------------------|
| 14044143 | hose & small phase separator |
| 14044151 | hose & medium phase separator |

CryoBio Accessory Catalog

Union Coupler

1110862

- used for connecting 2 LN2 transfer hoses together

Cryosleeve



9710899

100 per case

- keeps vials from falling from canes

Cabosil Block

11046952

small 2"

11026944

medium 3"

11046936

large 4"

- generally used in sets of 2
- used for CryoShipper Mini or to extend hold time of CryoShipper XC

CryoBio Accessory Catalog

Square Cabosil Insert



11046928

Secondary Shipping Container



10862106

- one part of the entire IATA package used for shipping bio-hazards

Screws for Secondary Shipping Container



2910221

CryoBio Accessory Catalog

Absorbant Paper for Secondary Shipping Container



10777446

100 sheets per bag

- each sheet absorbs 250ml
- used inside the secondary shipping container

CryoVials



10778810

1.2 ml vials

10778828

2 ml vials

10778836

3 ml vials

10777438

4 ml vials

- IATA approved shipping vials

Protective Shipping Containers

CryoBio Accessory Catalog

Protective Shipping Containers



| | |
|----------|---|
| 9719449 | SC 4/2V, SC 4/3V, SC 2/1V MiniMoover |
| 9722149 | XC 22/5 |
| 11925105 | SC 11/7, Doble 11 |
| 10537506 | CryoShipper, XC 20/3V |
| 11912460 | |
| 13027544 | CryoSystem 2000 w/ flat bottom XL thick foam 26"OD |
| 13027552 | XC 47/11, Doble 47 XL thick foam 26"OD |
| 11930861 | Doble 18 |
| 10741726 | CryoShipper extra capacity |

Roller Base for Aluminium Dewars



| | |
|----------|---|
| 10937391 | 14"ID 20 L units, Tech 35, Tech 50, SC 22/5 |
| 10937403 | 18.2"ID Lab 30, Lab 50, SC33/26, XC 34/18, MVE 230 |
| 10937411 | 20"ID XC47/11, MVE 450 |
| 10937420 | 22"ID XLC 360 HE, MVE 510, Cryo2000, Cryo4000 |
| 10937438 | 26.2"ID Cryo6000 |

Corks and Covers

CryoBio Accessory Catalog

Measuring Stick for Aluminum Dewars

5613469

Centimeters only (50 cm)
(24 inches)

Freezer Measuring Stick

10883361

1 meter in length
(39 inches)

Straw Canisters

CryoBio Accessory Catalog

Center Canister



| | |
|----------|-------------|
| 10654471 | CryoMoover |
| 10623201 | SC 20/20 |
| 10623181 | XC 47/11-6 |
| 10682181 | XC 22/5 |
| 10623199 | XC 34/18 |
| 10682181 | XC 20/3V |
| 10682199 | XC 32/8 |
| 11699347 | XC 47/11-10 |

C-1 Cane

9717009 1.0 ml 6/cane

- dimple but no dividers
- can hold goblets carrying straws

C-2 Cane



9717029 1.2 ml 6/cane

- no dividers or dimples
- can hold any length .5"OD vials or goblets carrying straws

CryoBio Accessory Catalog

C-7 Cane



9717089 2.0 ml 5/cane
- with dividers

Divider for Canes in Canister



10873286 XC 47/11-6 4 places
10873278 XC 47/11-6 6 places

Cabosil Canister Insert



10736409 1.50"OD
10736425 2.812"OD
10736441 3.968"OD
10736417 1.843"OD
10736433 3.375"OD
11501935 3.375"OD 6" length
11027719 3.187"OD

CryoBio Accessory Catalog

Plastic Mini Boxes




9719639
9719629

Holds 5x5 of 3" vial
Holds 5x5 of 2" vial

- holds internally or externally threaded vials

Large Plastic Vial Boxes



9179649

Holds 100 internally threaded vials 2" high

9719659

Holds 100 internally threaded vials 3" high

9719649CLR

Holds 100 Vials 2" with clear cover

All are 10"x10"

Large Cardboard Vial Boxes



9714099

Used for 81 externally threaded vials 2" high

9714109

used for 81 externally threaded vials 3" high

all are 9"x9"

CryoBio Accessory Catalog

Square Canister Rack



| | | |
|--|--------|-------------|
| 9719670 - XC 47/11-6 | Square | 3" Vial Box |
| 11554720 - 5 tier - will hold Nalgene vial boxes | Square | 2" Vial Box |
| 9719679 - CryoSystem 750 | Square | 3" Vial Box |

SUC Box



| | |
|----------|----------------------|
| 9710491 | |
| 11392644 | Rack for 2 SUC boxes |
| 11395896 | Rack for 8 SUC boxes |

SUC Marker Cane



10784574

CryoBio Accessory Catalog

SUC Cane Marker

| | |
|----------|--------|
| 11886679 | Silver |
| 11886687 | White |
| 11886695 | Brown |
| 11886716 | Pink |
| 11886724 | Blue |
| 11886732 | Green |
| 11886741 | Yellow |
| 11886759 | Orange |
| 11886767 | Red |

TEC 3000

13319459S

MVE Datalogger



| | |
|----------|---|
| 13820833 | Datalogger Only |
| 14075871 | Datalogger Kit |
| | -includes datalogger, software and comm cable |
| 14070202 | Cryoshipper, XC, and Mini |
| 14070288 | CryoMoover |
| 14070229 | SC 4/2V |
| 14070245 | SC 4/3V |
| 14070261 | XC 20/3V |
| 13792799 | Cryoshipper (c/c w/ fleece) |
| 13820884 | Replacement temp probe |
| 13822628 | Replacement mounting bracket |

CryoBio Accessory Catalog

Planer Datalogger



| | |
|----------|---------------------------|
| 11554324 | Datalogger only |
| 11371181 | Cryoshipper, XC, and Mini |
| 11554383 | CryoMoover |
| 11554341 | SC 4/2V |
| 11554367 | SC 4/3V |
| 11554375 | XC 20/3V |
| 11878353 | Replacement Planer probe |

Bat-1B Alarm



Liquid Level Alarm



| | |
|----------|-------------|
| 10769489 | Therm-o-let |
|----------|-------------|

CryoBio Accessory Catalog

1st Generation Conversion Kit



10740221 Converter Kit

includes:

- RJ45 cable assembly -8 conductor cable assembly with RJ45 plugs on each end
- DB25P to DB9S adaptor
- RJ45 to DB25S adaptor
- RS485 to RS232 converter
- Power supply and cord
- Daisy chain kit - RJ45 cable and splitter

Individual Part Numbers available upon request

2nd Generation Conversion Kit

13376947 USB Converter Kit

USB to RS-485

includes:

- converter
- RJ45 cable

Auto LN2 Tank Switcher

13934911

US tank switcher

11207040

Europe only tank switcher

Automatically switches liquid cylinders when supplying liquid nitrogen freezers. New high flow design to supply up to 8 freezers.

(does not include transfer hose)

CryoBio Accessory Catalog

Tee Valve Assembly



10784443

Used for connecting multiple supply tanks or multiple freezers

Vapor Platform



10544837

MVE 205

10544845

XLC 440

10544853

MVE 510

11509726

MVE 616

10544888

MVE 1223

10544888

MVE 1426

10544909

MVE 1839

Power Conditioner

11749885

assembled with mounting bracket

11807011

assembled without mounting bracket

Used for ensuring correct and constant 120V going to power supply without interference from power surges or lower voltage levels.

CryoBio Accessory Catalog

Freezer Step Assembly

| | |
|----------|--|
| 9721459 | foot step assembly 4 legged 1 step stool |
| 11507721 | free standing step 4 legged 1 step ss stool |
| 11662421 | single step mounts to 1500 series |
| 11822499 | two step mounts to 1500 series |
| 13082551 | single step mounts to 800 series |

Temperature Sleeves

| | |
|----------|----------|
| 11062055 | MVE 510 |
| 11937317 | MVE 616 |
| 10742518 | MVE 1223 |
| 11934044 | MVE 1426 |

Apron

10464394

CryoBio Accessory Catalog

Cryo-Gloves



Medium Length

| | |
|---------|------------------|
| 9717119 | medium size |
| 9717129 | large size |
| 9717139 | extra large size |

Elbow Length

| | |
|---------|------------------|
| 9717149 | medium size |
| 9717159 | large size |
| 9717169 | extra large size |

Cool Reach / Cryo Claw

| | |
|----------|----------------------|
| 13051579 | 40" reach |
| 13051608 | Replacement cups-set |

Epoxy Kit



| | |
|---------|--------------------------|
| 0010143 | 9 oz. resin and hardener |
|---------|--------------------------|

used for attaching canister ring to dewar or cork to cover