

NEW! INCUBATOR-GENIE™

BENCHTOP SHAKING/ROTATING INCUBATOR!



- Incubate, Rotate, Rock, and Shake all in ONE compact unit.
- Magnetic stainless steel platform may be easily set to rock or rotate almost any container.
- Precise temperature control and uniformity - range from 28°C to 75°C.
- Ideal for Hybridization.
- Chamber air is gently and continuously re-circulated ensuring temperature consistency and uniformity.
- Shatterproof window allows for inspection of samples without disturbing chamber environment.
- A variety of accessories are available to provide maximum application versatility.

Ordering Information

Catalog No.	Description
SI-1400	Incubator-Genie 120V
SI-1401	Incubator-Genie 230V, No Plug
SI-1402	Incubator-Genie 230V, European Plug
SI-1403	Incubator-Genie 230V, British Plug
SI-1404	Incubator-Genie 230V, Swiss Plug
SI-1406	Incubator-Genie 230V, Australian Plug
SI-1405	Incubator-Genie 100V

Specifications

Supplied	1 Clip Plate each for 10-13 mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf
Temperature	28°C - 75°C Accuracy: $\pm 0.2^\circ\text{C}$, Uniformity: $\pm 0.5^\circ\text{C}$ Resolution: 1.0°C , Stability: $\pm 0.1^\circ\text{C}$
Rotating	1-35 RPM
Rocking	2-70 cycles per minute
Platform Capacity	4.5 Kg (10 lb)
Chamber Volume	24 L (0.86 ft ³)
Weight	19.5 Kg (43lbs.)
Dimensions (D x W x H) (Overall)	400 x 560 x 370mm 15.75 x 22 x 14.5in.
Dimensions (D x W x H) (Chamber)	254 x 362 x 260mm 10 x 14.25 x 10.25in
Alarms	Visible and Audible for Temperature and Time



U.S. Pat. 5,918,979



Home of the Multi-Task Genie™ Family
70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: (888) 850-6208
Customer Service: custsvic@scientificindustries.com
www.scientificindustries.com

Incubator-Genie™ Accessories

Hold Virtually Any Container!



Accessories supplied with the Incubator-Genie




Fig.	Cat. No.	Description
1	SI-1123	1 Clip Plate for 12 each 10-13mm tubes
2	SI-1124	1 Clip Plate for 6 each 15-17mm tubes
3	SI-1125	1 Clip Plate for 3 each 28-30mm tubes
4	SI-1120	6 Clip Plates for 12 each 10-13mm tubes
5	SI-1121	6 Clip Plates for 6 each 15-17mm tubes
6	SI-1122	6 Clip Plates for 3 each 28-30mm tubes
7	SI-1126	Bag Mounting Strips (pack of 16)
8	SI-1127	Expansion Kit - two stainless steel trays, 24 bag mounting strips
9	SI-1130	2 Clip Plates for Hybridization Tubes
10	SI-1131	1 Stackable Shelf
11	SI-1134	2 Universal Clip Plates
12	SI-1135	1 Magnetic Covered 100 Microtube Box
13	SI-1136	1 Clip Plate for 100ml Volumetric Flasks
14	0A-1200-001	Non-Skid Tray
15	SI-1250	3-D Orbital Shaker Attachment

Plasticware, glassware and other consumable products shown are not included.



Scientific Industries, Inc. Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: (888) 850-6208
 Customer Service: custsvc@scientificindustries.com
 www.scientificindustries.com