




ABDOS CENTRIFUGE TUBES

Experience the Consistency & Purity with Abdos Centrifuge Tubes



Confirms US FDA 21 CFR, free from natural rubber, heavy metal & inhibitory plasticizers.

Even wall thickness to withstand maximum a force, 20,000 RCF.

Sterile tube have been validated according to ENISO 11137-2: SAL 10-6.

Precise lid sealing & perfect grip on screw for secure control on the cap in wet conditions.

Highly Transparent - 5% less haze percentage than a competitor.

Suitable for temperature of -80° to $+121^{\circ}$ C

Robust polymer selection - tubes do not leech & migrate unwanted substances into your sample.

Larger white frosted portion for easy writing.

Unique printing technology - chemical resistance, scratch proof for accurate & legible volume measurement.

Hydrodynamic Conical bottom for perfect peleting of sediments.

The Caps



Flat thread confirm vacuum tight closure - leak proof and prevent from override of cap.


Unique plug seal - ensuring perfect sealing & make it leak-proof.



P10430	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml Bulk, Blue Cap	500	25 Pcs / Bag
P10401	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Blue Cap	500	25 Pcs / Bag
P10441	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Red Cap	500	25 Pcs / Bag
P10402	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Sterile, Blue Cap	500	25 Pcs / Bag
P10442	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 15 ml, Sterile, Red Cap	500	25 Pcs / Bag
P10431	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Bulk, Blue cap	500	25 Pcs / Bag
P10403	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Blue Cap	500	25 Pcs / Bag

P10443	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Red Cap	500	25 Pcs / Bag
P10404	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Sterile, Blue Cap	500	25 Pcs / Bag
P10444	Abdos Centrifuge Tubes Conical Bottom, PP /HOPE, 50 ml, Sterile, Red Cap	500	25 Pcs / Bag
P10405	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 15 ml, Amber, Blue Cap	500	25 Pcs / Bag
P10406	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 15 ML, Sterile Amber, Blue Cap	500	25 Pcs / Bag
P10407	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 50 ml, Amber, Blue Cap	500	25 Pcs / Bag
P10408	Abdos Centrifuge Tubes Conical Bottom, PP/ Amber, 50 ml, Sterile Amber, Blue Cap	500	25 Pcs / Bag
P10432	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Bulk, Blue Cap	500	25 Pcs / Bag
P10423	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Blue Cap	500	25 Pcs / Bag
P10449	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Red Cap	500	25 Pcs / Bag
P10424	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Sterile, Blue Cap	500	25 Pcs / Bag
P10450	Abdos Self Standing Centrifuge Tubes, PP/HOPE, 50 ml, Sterile, Red Cap	500	25 Pcs / Bag
P10419	Abdos Centrifuge Tubes Conical Bottom, PS /HOPE, 15 ml, Sterile, Blue Cap	500	25 Pcs / Bag
P10421	Abdos Centrifuge Tubes Conical Bottom, PS /HOPE, 15 ml, Racked Sterile, Blue Cap	750	25 Pcs / Bag

P10425	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 15 ml, Sterile, Blue Cap	500	1 Ps / Bag
P10451	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 15 ml, Sterile, Red Cap	500	1 Ps / Bag
P10426	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Blue Cap	500	1 Ps / Bag
P10452	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Red Cap	500	1 Ps / Bag
P10427	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Self Standing, Blue Cap	450	1 Ps / Bag
P10453	Abdos Individually Wrapped Centrifuge Tubes, PP /HOPE, 50 ml, Sterile, Self Standing, Red Cap	450	1 Ps / Bag



Centrifuge Tube 50 ml
(Conical Bottom)



Centrifuge Tube 50 ml
(Self stand)



Centrifuge Tube 15 ml
(Conical bottom)